EasyDisc* R2A Test Kit

06-0013525-01

IDEXX

For Technical Support, please call:

North/South America: 1 207 556 4496/1 800 321 0207

Europe: 00800 4339 9111

UK: +44 (0) 1638 676800

China: +86-21-61279528

Japan: 03 5301 6800

Australia: 1300 443 399

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

idexx.com/water

EasyDisc® R2A Test Kit

Introduction

The EasyDisc R2A method is used for the quantification of culturable, waterborne bacteria. Culturable is defined as any aerobic bacteria able to form colonies on solid media. It is based on colorimetric technology which detects viable aerobic bacteria by testing for growth and for the presence of key enzymes known to be present in these organisms. It uses a substrate that produces a blue color reaction when metabolized by most waterborne bacteria. The sample is added directly to the EasyDisc R2A plate, incubated, and then examined for blue colonies. The colony values generated by EasyDisc R2A method correlate with the Pour Plate method using Reasoners 2 Agar (R2A) incubated at 20-28°C for 5 to 7 days as described in *Standard Methods for the Examination of Water and Wastewater: SM9215 Heterotrophic Plate Count (2017)*¹.

Contents

25 or 100 sterile EasyDisc R2A plates in plastic sleeves.

Storage

Store at 2-25°C away from light and humidity. Expiration date is located on the outer kit label and Quality Control Certificate.

Test Procedure

1. Remove white lid and add 1 ± 0.1 mL water sample directly to the clear EasyDisc R2A plate.
2. Gently swirl *immediately* after sample addition to coat the plate entirely and reattach white lid.
3. Incubate undisturbed at room temperature for at least 20 minutes to allow for the media to set. Transfer plate to incubator within 1 hour of sample addition.
4. Incubate the plate, white lid side up, at 20-28°C for 5 to 7 days.
5. Examine the plates for colony growth as soon as they are removed from the incubators. Reject any plate with confluent growth.
6. Sum all colonies for the result. The result is expressed as colony forming units (CFU) per milliliter.

NOTE: Most microorganisms will produce a blue color on EasyDisc, but some may produce natural pigments with colors other than blue. All colonies should be summed for the final result.

Procedural Notes

1. This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
2. Chlorinated samples should be treated with sodium thiosulfate prior to testing, preferably upon sample collection.
3. Printed grid is designed to aid in counting colonies.
4. Follow aseptic technique. Dispose of sample and media in accordance with Good Laboratory Practice.
5. Samples may be diluted before adding to the media provided the final volume applied is 1 ± 0.1 mL. The recommended sterile diluents are dechlorinated water and deionized water.
6. Adjust the CFU/mL result to reflect dilutions. For example, if a 10-fold dilution is prepared by adding 0.1 mL of sample into 0.9 mL of diluent, multiply result by 10 to convert to CFU/1 mL.
7. If there are more than 300 colonies on the plates inoculated with the highest dilutions used, express the results as >300 CFU/mL or approximate only.

Quality Control Procedures

1. The following quality control procedures are recommended for each lot of EasyDisc R2A:
 - a. Positive Control option 1: IDEXX-QC HPC/TVC²: *Enterococcus faecalis*
 - b. Positive Control option 2: Fill a sterile vessel with 100 mL of sterile, non-buffered, oxidant-free water and inoculate with a sterile loop of any of the following strains:

Organism	WDCM#	ATCC#
<i>Enterococcus faecalis</i>	00087	29212
<i>Escherichia coli</i>	00090	11775

- c. Negative Control/Blank: Use 1 ml of sterile diluent.
2. Follow Test Procedure Steps 1-6.
 3. Negative Control/Blank test should not contain any colonies.

NOTE: IDEXX internal quality control testing is performed in accordance with ISO 11133:2014³. Quality Control Certificates are available at idexx.com/water.

The EasyDisc R2A Test is intended to be used only for water quality research and analysis, by technically qualified individuals or under their supervision.

1. *Standard Methods for the Examination of Water and Wastewater: SM9215 Heterotrophic Plate Count (2017)*.

2. IDEXX-QC HPC/TVC, IDEXX Catalog #UN3373-WQC-HPC

3. ISO International Standardization Organisation. *Microbiology of food, animal feed and water - Preparation, production, storage and performance testing of culture media*. ISO 11133:2014.

*EasyDisc is a registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. Patent information: idexx.com/patents.

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¹ Standard Methods for the Examination of Water and Wastewater: SM9215 Heterotrophic Plate Count (2017).

² IDEXX-QC HPC/TVC, IDEXX Catalog #UN3373-WOC-HPC

³ ISO International Standardization Organisation. *Microbiology of food, animal feed and water - Preparation, production, storage and performance testing of culture media*. ISO 11133:2014.

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1. Standard Methods for the Examination of Water and Wastewater: SM9215 Heterotrophic Plate Count (2017).

2. IDEXX-QC HPC/TVC, IDEXX Catalog # UN3373-WQC-HPC

3. ISO International Standardization Organisation. Microbiology of food, animal feed and water - Preparation, production, storage and performance testing of culture media. ISO 11133:2014.

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© 2021 IDEXX Laboratories, Inc. 無断複写・転載を禁じます。