
Enterolert^{*} 250

06-26891-04

IDEXX

For Technical Support, please call:

North/South America: 1 207 556 4496/1 800 321 0207

Europe: 00800 4339 9111

UK: +44 (0) 1638 676800

China: +86-21-61279528

Japan: +03 5301 6800

Australia: 1300 443 399

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

idexx.com/water

Enterolert* 250 Test Kit

Introduction

Enterolert* 250 is designed to detect *Enterococcus* species in bottled water samples. It is based on IDEXX's proprietary Defined Substrate Technology* (DST*). Enterolert 250 provides confirmed presence/absence results in 24 hours. Enterolert 250 uses Ortho-Nitrophenyl- β -D-Glucoside as a nutrient indicator and incorporates a specifically designed blue background colour in its formulation. When the substrate is metabolized by enterococci, the sample turns from blue to green to indicate a positive result. Any change from the original colour to green is considered a positive result. No ultraviolet light source is required. Enterolert 250 detects enterococci at 1 cfu in 250 mL bottled water samples in 24 hours.

Storage

Store at 2–25°C away from direct light and humidity.


Presence/Absence (P/A) Procedure for Noncarbonated Water Samples

1. Add contents of one snap pack to a 250 mL sample in a sterile, transparent vessel.
2. Cap vessel and shake until contents are dissolved.
3. Incubate at $41 \pm 1^\circ\text{C}$ for 24 hours; results are valid for up to 28 hours of incubation.
4. Read results according to Result Interpretation table, below.

Presence/Absence (P/A) Procedure for Carbonated Water Samples

1. De-gas sample by vigorously shaking sample in its original container until no gas can be seen or heard (loosen cap to allow gas to escape after shaking).
2. Add contents of one snap pack to a 250 mL sample in a sterile, transparent vessel.
3. Shake vessel vigorously for 30–60 seconds to release carbonation (loosen cap to allow gas to escape after shaking).
4. Secure cap and incubate at $41 \pm 1^\circ\text{C}$ for 24 hours; results are valid for up to 28 hours of incubation.
5. Read results according to Result Interpretation table, below.

Result Interpretation

Appearance	Result
Blue	Negative for enterococci
Green	Positive for enterococci

- Enterolert 250 results are definitive at 24–28 hours. In addition, positives for enterococci observed before 24 hours and negatives observed after 28 hours are also valid.

Procedural Notes

- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a negative control can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Enterolert 250 is a primary water test. Enterolert 250 performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- Enterolert 250 has not been validated for use with flavored bottled water.
- The presence of a high mineral content (especially magnesium and/or calcium) can cause Enterolert 250 reagent to become cloudy but does not affect test result outcome.
- Aseptic technique should always be followed when using Enterolert 250. Dispose of in accordance with Good Laboratory Practice.

Quality Control Procedures

1. One of the following quality control procedures is suggested for each lot of Enterolert 250:
 - A. IDEXX-QC Enterococci DW: *Enterococcus faecalis*, *Serratia marcescens*
 - B.
 - i. For each of the American Type Culture Collection² (ATCC) bacterial strains (*Enterococcus faecalis* 29212, *Aerococcus viridans* 11563), streak the culture onto labeled TSA or blood agar plates and incubate at $35 \pm 2^\circ\text{C}$ for 18–24 hours.
 - ii. For positive controls, prepare proper dilutions to reach a concentration of approximately 10–50 CFU/250 mL.
 - iii. For non-target strains, prepare proper dilutions to a suggested inoculation up to 10^6 CFU/250 mL.
2. Follow the P/A Procedure above.
3. Results should match the Result Interpretation table above.

1. IDEXX-QC Enterococci DW. IDEXX Catalog #UN3373-WQC ENT DW
2. American Type Culture Collection 1-800-638-6597 www.atcc.org

*Defined Substrate Technology, DST and Enterolert are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.
Patent information: idexx.com/patents.

© 2019 IDEXX Laboratories, Inc. All rights reserved.

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© 2019 IDEXX Laboratories, Inc. Wszelkie prawa zastrzeżone.